

Product Datasheet

Human Proteolipid Protein 1, Myelin (PLP1) ELISA Kit (orb1207836)

Description	Human Proteolipid Protein 1, Myelin (PLP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P60201
Sensitivity	0.053ng/mL
Entrez	5354
Expiration Date	Please enquire.