

Product Datasheet

Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2a) ELISA Kit (orb1207850)

Description	Human Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide (PIK3C2a) ELISA Kit
Reactivity	Human
Range	3.12-200U/L
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	O00443
Sensitivity	1.6U/L
Entrez	5286
Expiration Date	Please enquire.