

## Product Datasheet

# Human Poly ADP Ribose Polymerase (PARP) ELISA Kit (orb1207872)

<b>Description</b>	Human Poly ADP Ribose Polymerase (PARP) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.781-50ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P09874</a>
<b>Sensitivity</b>	0.25ng/mL
<b>Entrez</b>	<a href="#">142</a>
<b>Expiration Date</b>	Please enquire.