

Product Datasheet

Human Meprin A Alpha (MEP1a) ELISA Kit (orb1207935)

Description	Human Meprin A Alpha (MEP1a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q16819
Sensitivity	0.057ng/mL
Entrez	4224
Expiration Date	Please enquire.