

Product Datasheet

Human Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit (orb1207937)

Description	Human Monocyte Chemotactic Protein 2 (MCP2) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P80075
Sensitivity	8.6pg/mL
Entrez	6355
Expiration Date	Please enquire.