

Product Datasheet

Human Inhibin Beta E (INHbE) ELISA Kit (orb1207977)

Description	Human Inhibin Beta E (INHbE) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P58166
Sensitivity	6.4pg/mL
Entrez	83729
Expiration Date	Please enquire.