

Product Datasheet

Human Inducible T-Cell Co Stimulator Ligand (ICOSLG) ELISA Kit (orb1207981)

Description	Human Inducible T-Cell Co Stimulator Ligand (ICOSLG) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O75144
Sensitivity	0.057ng/mL
Entrez	23308
Expiration Date	Please enquire.