

## Product Datasheet

# Human High Mobility Group Protein 1 (HMG1) ELISA Kit (orb1208001)

<b>Description</b>	Human High Mobility Group Protein 1 (HMG1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5-4000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<b>P09429</b>
<b>Sensitivity</b>	23.7pg/mL
<b>Entrez</b>	<b>3146</b>
<b>Expiration Date</b>	Please enquire.