

Product Datasheet

Human Glutathione S Transferase Mu 1 (GSTm1) ELISA Kit (orb1208015)

| | |
|--------------------------|---|
| Description | Human Glutathione S Transferase Mu 1 (GSTm1) ELISA Kit |
| Reactivity | Human |
| Range | 1.56-100ng/mL |
| Note | For research use only |
| Application notes | sandwich |
| Sample Types | serum, plasma, tissue homogenates or other biological fluids. |
| Uniprot ID | P09488 |
| Sensitivity | 0.52ng/mL |
| Entrez | 2944 |
| Expiration Date | Please enquire. |