

## Product Datasheet

# Human CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit (orb1208119)

<b>Description</b>	Human CCAAT/Enhancer Binding Protein Beta (CEBPb) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<b>P17676</b>
<b>Sensitivity</b>	0.061ng/mL
<b>Entrez</b>	<b>1051</b>
<b>Expiration Date</b>	Please enquire.