

## Product Datasheet

# Human Chemokine C-C-Motif Receptor 2 (CCR2) ELISA Kit (orb1208133)

<b>Description</b>	Human Chemokine C-C-Motif Receptor 2 (CCR2) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P41597</a>
<b>Sensitivity</b>	0.061ng/mL
<b>Entrez</b>	<a href="#">729230</a>
<b>Expiration Date</b>	Please enquire.