

## Product Datasheet

# Human Chemokine C-C-Motif Ligand 3 Like Protein 1 (CCL3L1) ELISA Kit (orb1208136)

<b>Description</b>	Human Chemokine C-C-Motif Ligand 3 Like Protein 1 (CCL3L1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.1-5000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P16619</b>
<b>Sensitivity</b>	32pg/mL
<b>Entrez</b>	<b>6349</b>
<b>Expiration Date</b>	Please enquire.