

Product Datasheet

Human Aquaporin 1, Colton Blood Group (AQP1) ELISA Kit (orb1208172)

Description	Human Aquaporin 1, Colton Blood Group (AQP1) ELISA Kit
Reactivity	Human
Range	0.25-16ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P29972
Sensitivity	0.09ng/mL
Entrez	358
Expiration Date	Please enquire.