

Product Datasheet

Human Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit (orb1208181)

Description	Human Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	6.43pg/mL
Expiration Date	Please enquire.