

Product Datasheet

Rat Interleukin 33 (IL33) ELISA Kit (orb1208224)

Description	Rat Interleukin 33 (IL33) ELISA Kit
Reactivity	Rat
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q66H70
Sensitivity	27pg/mL
Entrez	361749
Expiration Date	Please enquire.