

Product Datasheet

Canine Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit (orb1208252)

Description	Canine Interferon Gamma Induced Protein 10kDa (IP10) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q5KSV9
Sensitivity	6.8pg/mL
Entrez	478432
Expiration Date	Please enquire.