

Product Datasheet

Canine Interleukin 5 (IL5) ELISA Kit (orb1208257)

Description	Canine Interleukin 5 (IL5) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q95J76
Sensitivity	5.4pg/mL
Entrez	403790
Expiration Date	Please enquire.