

Product Datasheet

Canine Interleukin 13 (IL13) ELISA Kit (orb1208275)

Description	Canine Interleukin 13 (IL13) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NOW9
Sensitivity	5.5pg/mL
Entrez	442990
Expiration Date	Please enquire.