

Product Datasheet

Canine Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) ELISA Kit (orb1208286)

Description	Canine Colony Stimulating Factor 2, Granulocyte Macrophage (GMCSF) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48749
Sensitivity	7.2pg/mL
Entrez	403923
Expiration Date	Please enquire.