

## Product Datasheet

# Guinea pig Endothelin 1 (EDN1) ELISA Kit (orb1208295)

<b>Description</b>	Guinea pig Endothelin 1 (EDN1) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	7.81-500pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>P97740</b>
<b>Sensitivity</b>	2.82pg/mL
<b>Entrez</b>	<b>100726197</b>
<b>Expiration Date</b>	Please enquire.