

Product Datasheet

Canine Endothelin 1 (EDN1) ELISA Kit (orb1208296)

Description	Canine Endothelin 1 (EDN1) ELISA Kit
Reactivity	Canine
Range	7.81-500pg/mL
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13206
Sensitivity	2.33pg/mL
Entrez	403424
Expiration Date	Please enquire.