

## Product Datasheet

# General Tetrahydrobiopterin (THB) ELISA Kit (orb1208310)

<b>Description</b>	General Tetrahydrobiopterin (THB) ELISA Kit
<b>Reactivity</b>	Other
<b>Range</b>	0.156-10ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates or other biological fluids.
<b>Sensitivity</b>	0.056ng/mL
<b>Expiration Date</b>	Please enquire.