

## Product Datasheet

# General Prostacyclin (PGI<sub>2</sub>) ELISA Kit (orb1208314)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | General Prostacyclin (PGI <sub>2</sub> )<br>ELISA Kit |
| <b>Reactivity</b>        | Other   |
| <b>Range</b>             | 78.1-5000pg/mL  |
| <b>Note</b>              | For research use only                                 |
| <b>Application notes</b> | competitive inhibition                                |
| <b>Sample Types</b>      | serum, plasma or other<br>biological fluids.          |
| <b>Sensitivity</b>       | 33.5pg/mL   |
| <b>Expiration Date</b>   | Please enquire.                                       |