

## Product Datasheet

# General Prostaglandin F2 Alpha (PGF2a) ELISA Kit (orb1208317)

<b>Description</b>	General Prostaglandin F2 Alpha (PGF2a) ELISA Kit
<b>Reactivity</b>	Other
<b>Range</b>	31.2-2000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	11.67pg/mL
<b>Expiration Date</b>	Please enquire.