

## Product Datasheet

# General 25-Hydroxyvitamin D3 (HVD3) ELISA Kit (orb1208321)

<b>Description</b>	General 25-Hydroxyvitamin D3 (HVD3) ELISA Kit
<b>Reactivity</b>	Other
<b>Range</b>	3.12-200ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Sensitivity</b>	0.83ng/mL
<b>Expiration Date</b>	Please enquire.