

Product Datasheet

Human Growth Factor Receptor Bound Protein 10 (Grb10) ELISA Kit (orb1208408)

Description	Human Growth Factor Receptor Bound Protein 10 (Grb10) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q13322
Sensitivity	0.102ng/mL
Entrez	2887
Expiration Date	Please enquire.