

Product Datasheet

Mouse Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) ELISA Kit (orb1208432)

Description	Mouse Tumor Necrosis Factor Ligand Superfamily, Member 13 (TNFSF13) ELISA Kit
Reactivity	Mouse
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9D777
Sensitivity	22.7pg/mL
Entrez	69583
Expiration Date	Please enquire.