

## Product Datasheet

# Mouse Interleukin 5 Receptor Alpha (IL5Ra) ELISA Kit (orb1208464)

<b>Description</b>	Mouse Interleukin 5 Receptor Alpha (IL5Ra) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	62.5-4000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">P21183</a>
<b>Sensitivity</b>	28.5pg/mL
<b>Entrez</b>	<a href="#">16192</a>
<b>Expiration Date</b>	Please enquire.