

Product Datasheet

Chicken Interferon Alpha (IFNa) ELISA Kit (orb1208521)

Description	Chicken Interferon Alpha (IFNa) ELISA Kit
Reactivity	Gallus
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P42165
Sensitivity	28pg/mL
Entrez	100857744
Expiration Date	Please enquire.