

Product Datasheet

Rat Endothelin 1 (EDN1) ELISA Kit (orb1208542)

Description	Rat Endothelin 1 (EDN1) ELISA Kit
Reactivity	Rat
Range	7.81-500pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P22388
Sensitivity	2.19pg/mL
Entrez	24323
Expiration Date	Please enquire.