

Product Datasheet

Rabbit Connective Tissue Growth Factor (CTGF) ELISA Kit (orb1208543)

Description	Rabbit Connective Tissue Growth Factor (CTGF) ELISA Kit
Reactivity	Rabbit
Range	46.8-3000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q2WG38
Sensitivity	18.54pg/mL
Entrez	100008638
Expiration Date	Please enquire.