

Product Datasheet

Mouse BH3 Interacting Domain Death Agonist (Bid) ELISA Kit (orb1208584)

Description	Mouse BH3 Interacting Domain Death Agonist (Bid) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P70444
Sensitivity	0.112ng/mL
Entrez	12122
Expiration Date	Please enquire.