

Product Datasheet

Porcine Interleukin 8 (IL8) ELISA Kit (orb1208590)

Description	Porcine Interleukin 8 (IL8) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P26894
Sensitivity	12.3pg/mL
Entrez	396880
Expiration Date	Please enquire.