

Product Datasheet

Human Interleukin 18 Receptor 1 (IL18R1) ELISA Kit (orb1208629)

Description	Human Interleukin 18 Receptor 1 (IL18R1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q13478
Sensitivity	0.058ng/mL
Entrez	8809
Expiration Date	Please enquire.