

Product Datasheet

Human Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit (orb1208633)

Description	Human Interleukin 13 Receptor Alpha 2 (IL13Ra2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q14627
Sensitivity	0.129ng/mL
Entrez	3598
Expiration Date	Please enquire.