

Product Datasheet

Human Phosphatase And Tensin Homolog (PTEN) ELISA Kit (orb1208661)

Description	Human Phosphatase And Tensin Homolog (PTEN) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P60484
Sensitivity	0.136ng/mL
Entrez	5728
Expiration Date	Please enquire.