

Product Datasheet

Human Heat Shock 70kDa Protein 8 (HSPA8) ELISA Kit (orb1208672)

Description	Human Heat Shock 70kDa Protein 8 (HSPA8) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P11142
Sensitivity	0.134ng/mL
Entrez	3312
Expiration Date	Please enquire.