

Product Datasheet

Mouse Transforming Growth Factor Beta Receptor I (TGFbR1) ELISA Kit (orb1208708)

Description	Mouse Transforming Growth Factor Beta Receptor I (TGFbR1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q64729
Sensitivity	0.069ng/mL
Entrez	21812
Expiration Date	Please enquire.