

Product Datasheet

Human Interleukin 1 Receptor Type I (IL1R1) ELISA Kit (orb1208793)

Description	Human Interleukin 1 Receptor Type I (IL1R1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P14778
Sensitivity	0.049ng/mL
Entrez	3554
Expiration Date	Please enquire.