

## Product Datasheet

# General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit (orb1208825)

<b>Description</b>	General Cyclic Adenosine Monophosphate (cAMP) ELISA Kit
<b>Reactivity</b>	Other
<b>Range</b>	0.312-20ng/mL
<b>Note</b>	For research use only
<b>Application notes</b>	competitive inhibition
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	0.118ng/mL
<b>Expiration Date</b>	Please enquire.