

## Product Datasheet

# Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit (orb1208829)

<b>Description</b>	Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	46.8-3000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q7TNL0</a>
<b>Sensitivity</b>	16.13pg/mL
<b>Entrez</b>	<a href="#">305236</a>
<b>Expiration Date</b>	Please enquire.