

Product Datasheet

Mouse Alpha-2-Macroglobulin (a2M) ELISA Kit (orb1208839)

Description	Mouse Alpha-2-Macroglobulin (a2M) ELISA Kit
Reactivity	Mouse
Range	9.37-600ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, urine or other biological fluids.
Uniprot ID	Q6GQT1
Sensitivity	3.7ng/mL
Entrez	232345
Expiration Date	Please enquire.