

Product Datasheet

Human Macrophage Inflammatory Protein 4 Alpha (MIP4a) ELISA Kit (orb1208862)

Description	Human Macrophage Inflammatory Protein 4 Alpha (MIP4a) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9Y258
Sensitivity	6.1pg/mL
Entrez	10344
Expiration Date	Please enquire.