

Product Datasheet

Guinea pig Tumor Necrosis Factor Alpha (TNFa) ELISA Kit (orb1208877)

Description	Guinea pig Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P51435
Sensitivity	5.5pg/mL
Entrez	100135630
Expiration Date	Please enquire.