

## Product Datasheet

# Guinea pig Interleukin 1 Beta (IL1b) ELISA Kit (orb1208884)

<b>Description</b>	Guinea pig Interleukin 1 Beta (IL1b) ELISA Kit
<b>Reactivity</b>	Guinea pig
<b>Range</b>	78.1-5000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9WVG1</a>
<b>Sensitivity</b>	25pg/mL
<b>Entrez</b>	<a href="#">100135556</a>
<b>Expiration Date</b>	Please enquire.