

Product Datasheet

Canine Interleukin 1 Beta (IL1b) ELISA Kit (orb1208885)

Description	Canine Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q2VCE2
Sensitivity	6.4pg/mL
Entrez	403974
Expiration Date	Please enquire.