

## Product Datasheet

# Human Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) ELISA Kit (orb1208888)

<b>Description</b>	Human Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6-1000pg/mL
<b>Note</b>	For research use only
<b>Application notes</b>	sandwich
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Uniprot ID</b>	<b>Q99075</b>
<b>Sensitivity</b>	6.3pg/mL
<b>Entrez</b>	<b>1839</b>
<b>Expiration Date</b>	Please enquire.