

Product Datasheet

Mouse Macrophage Inflammatory Protein 3 Alpha (MIP3a) ELISA Kit (orb1208897)

Description	Mouse Macrophage Inflammatory Protein 3 Alpha (MIP3a) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	O89093
Sensitivity	27pg/mL
Entrez	20297
Expiration Date	Please enquire.