

Product Datasheet

Mouse Chemerin (CHEM) ELISA Kit (orb1208904)

Description	Mouse Chemerin (CHEM) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9DD06
Sensitivity	26pg/mL
Entrez	71660
Expiration Date	Please enquire.