

Product Datasheet

Rat Galectin 1 (GAL1) ELISA Kit (orb1208939)

Description	Rat Galectin 1 (GAL1) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11762
Sensitivity	0.246ng/mL
Entrez	56646
Expiration Date	Please enquire.